




Memorial Sloan Kettering
Cancer Center™

Molecular Imaging of Ovarian Cancer



<http://www.nature.com/nmeth/journal/v6/n6/images/nmeth0609-465-l2.jpg>


Robert Kimelman
Secaucus High School – Secaucus, NJ
Summer Student Intern
Lewis Lab Group Meeting
August 26, 2016




<http://www.genengnews.com/gen-articles/the-antibody-dilemma/5550/>

Background Information

- Symptoms: very subtle and non-specific
- Most fatal gynecologic malignancy
- Molecular signature of LI CAM antigen
- Sought to noninvasively perceive LI CAM in SKOV3 cell lines via PET/CT



<http://ovariancancerday.org/about-ovarian/>





<http://www.dallas-obgyn.com/gynecology/ovarian-cancer/>



Background Information

- Positron emission tomography is at forefront of cancer “theranostics”
- Detects gamma rays emitted by isotopes to visualize radioactivity
- PET radiopharmaceuticals consist of:
 - A molecular structure
 - A positron emitting radionuclide
- Zirconium-89 labeling is favorable for PET
- Use of desferrioxamine (DFO) is very effective



https://www.researchgate.net/figure/262763777_fig1_Schematic-overview-of-89Zr-labeled-antibody-using-DFO-as-chelator

<http://www.sciencedirect.com/science/article/pii/S0720048X09006913>



^{89}Zr -DFO-H3E7 **IgG4 Mutated**

24h






48h

72h






96h

120h


Mouse 1



Mouse 2



H3E7-4M 20h Biodistrib; on&



Conjugation/Labeling Conditions

H3E7-4M DFO Conjugation

- 500 ug antibody
- DFO-antibody molar ratio – 20:1
- pH 9, 37 °C for 60 minutes in 1 mL

⁸⁹ZrDFO-H3E7-4M Labeling


- 236 ug antibody, 650 uCi Zr
- pH 7.4, RT for 60 minutes
- **Crude RCY: 95% RCP: 99% SA: 2.75 uCi/ug**




Serum Stability




Immunoreactivity



FACS (cell line: IMR32-Luc)



Zirconium-89 Labeled H3E7-4M Internalization Assay



Conclusions and Future Work

- Minimal uptake of Zirconium-89 labeled H3E7-4M seen in kidneys
 - Significant tumor uptake of Zirconium-89 labeled H3E7-4M
 - SKOV3+ cells are actively internalizing Zirconium-89 labeled H3E7-4M
 - Zirconium-89 labeled H3E7-4M is about 80% immunoreactive
-
- Future Work:
 - Images with higher specific activity and biodistribution



Techniques I've Learned

- Cell Culture
- Pipetting
- PET/CT Imaging
- HPLC
- Western Blotting
- Antibody Conjugation
- Radiolabeling
- Biodistribution
- RARC Certified
- Tube Weighing
- Mouse Subcutaneous Xenografting
- Mouse Intravenous Tail Injections
- Immunoreactivity Cell Assay
- Internalization Cell Assay
- Zirconium-89 Preparation
- Orthotopic (Mouse) Surgery



Acknowledgements

- Brandon Nemieboka
- Dr. Jason Lewis
- The Lewis Lab

